


ADAPTERS

QA-1980 ADAPTER VACUUM, STRAIGHT

An outer joint at top with inner joint of the same size at bottom. Delivery tube is sealed in above hose connection.


<u>Joint Size</u>	<u>Code</u>
24/40	-110
29/42	-111



QA-2020 ADAPTER VACUUM, 105°

An outer joint at top on a 105° angle with inner joint of the same size at bottom. Delivery tube is sealed in above a side hose connection. Please specify tube length below joint.

<u>Joint Size</u>	<u>Code</u>
24/40	-110
29/42	-111



QA-2000 ADAPTER GAS INLET

A right angle inlet tube extends through an inner joint with side hose connection. Please specify tube length below joint.

<u>Joint Size</u>	<u>Code</u>
19/38	-110
24/40	-111
29/42	-112



REMEMBER

“Q” Glass has been working with customers to design and build custom glass apparatus and component parts since 1946. We can do the same for you. Please write or call us with your special requirements.

ARSENIC APPARATUS

QA-3000 ARSENIC APPARATUS

Consists of a 125ml erlenmeyer flask with a 24/40 outer joint. A scrubber tube with a 24/40 inner joint and 12/2 ball joint, an absorber tube with a 12/2 socket joint and a #12 clamp.

Code
-110

