

BULBS

QB-1180 BULB KJELDAHL

An outer joint at the top with inner joint below a round trap.


| Joint Size | Code |
|------------|------|
| 24/40 | -110 |
| 29/42 | -111 |



QB-1220 BULB KJELDAHL

Inner joints on the bottom and side arm are the same size. Bottom inner tube has opening for return of trapped liquid.


| Bulb Size | | |
|------------|-------|-------|
| Diameter/ | 24/40 | 29/42 |
| Length, MM | Code | Code |
| 45/100 | -110 | -212 |
| 55/125 | -111 | -213 |



QB-1200 BULB KJELDAHL

A drip joint below a round trap with curved inner tube on 75° side arm.

| Joint Size | Code |
|------------|------|
| 24/40 | -110 |
| 29/42 | -111 |



REMEMBER

“Q” Glass has been working with customers to design and build custom glass apparatus and component parts since 1946. We can do the same for you. Please write or call us with your special requirements.

BURETS

QB-1260 BURET GLASS STOPCOCK


Laboratory grade buret, graduated, and has a 2mm straight stopcock.

| Capacity, ML | Subdivisions, ML | Code |
|--------------|------------------|------|
| 10 | 0.05 | -110 |
| 25 | 0.1 | -111 |
| 50 | 0.1 | -112 |
| 100 | 0.2 | -113 |

QB-1280 BURET TEFLON STOPCOCK

Similar to QB-1260, but with a 2mm teflon stopcock.

| Capacity, ML | Subdivisions, ML | Code |
|--------------|------------------|------|
| 10 | 0.05 | -110 |
| 25 | 0.1 | -111 |
| 50 | 0.1 | -112 |
| 100 | 0.2 | -113 |



QB-1300 BURET GLASS STOPCOCK

Similar to QB-1260, but with the addition of a schellbach stripe.

| Capacity, ML | Subdivisions, ML | Code |
|--------------|------------------|------|
| 10 | 0.05 | -110 |
| 25 | 0.1 | -111 |
| 50 | 0.1 | -112 |
| 100 | 0.2 | -113 |

QB-1320 BURET TEFLON STOPCOCK

Similar to QB-1280, but with the addition of a schellbach stripe.

| Capacity, ML | Subdivisions, ML | Code |
|--------------|------------------|------|
| 10 | 0.05 | -110 |
| 25 | 0.1 | -111 |
| 50 | 0.1 | -112 |
| 100 | 0.2 | -113 |

