


CLAMPS

QC-2011 CLAMP HOLDER

To hold clamps to support rods.
Clamp will grip rods up to 19mm.


Code
-110



QC-2013 TRIPOD SUPPORT

Base is heavy cast iron for stability, chemical resistant black enamel. Rod and base are precisely threaded.


<u>Rod Length and Diameter, In.</u>	<u>Code</u>
18 x 5/16	-110
20 x 3/8	-111
24 x 7/16	-112
36 x 1/2	-113



QC-2012 CLAMP RING SUPPORT

2" split ring allows easy placement and removal of objects. Vinyl coated for protection.

Code
-110



QC-2014 TUBING CLAMP

Adjusts easily with a screwdriver. Made of stainless steel. Sold in packages of 10 each.

<u>Fits Tubing Sizes, OD</u>	<u>Code</u>
1/4" - 5/8"	-110
5/16" - 7/8"	-111




CLINICAL APPARATUS

QC-2015 WARBURG REACTION VESSEL

Approximately 15ml with center well and 7/15 solid stopper in side arm.
Center joint is a 17/20.


Code
-110



QC-2020 WARBURG REACTION VESSEL

Similar to QC-2015 but with 7/15 venting plug in side arm.

Code
-110



REMEMBER


"Q" Glass has been working with customers to design and build custom glass apparatus and component parts since 1946. We can do the same for you. Please write or call us with your special requirements.

CLINICAL APPARATUS

QC-2030 WARBURG REACTION VESSEL

Approximately 15ml without center well, 7/15 solid stopper in side arm. Center joint is a 17/20.


Code
-110



QC-2050 WARBURG MANOMETER

Twin capillaries graduated 0 to 30cm in 1mm divisions, easily read. Tail outlet stopcock on top and a 17/20 joint on side arm for attaching reaction vessels QC-2015 through QC-2040.


Code
-110



QC-2040 WARBURG REACTION VESSEL

Approximately 15ml with center well and 2 side arms. One side arm has a 7/15 solid stopper. The other side arm has a 7/15 venting pipe.

Code
-110




CONDENSERS

QC-2080 CONDENSER ALLIHN

The top of the condenser is a plain tube. The bottom has a drip tip.

Jacket Length	24/40	29/42
MM	<u>Code</u>	<u>Code</u>
200	-110	-216
250	-111	-217
300	-112	-218
400	-113	-219
500	-114	-220
600	-115	-221



QC-3000 CONDENSER ALLIHN

Similar to QC-2080, but with an outer joint at top and a drip tip of the same size at bottom.

Jacket Length	24/40	29/42
MM	<u>Code</u>	<u>Code</u>
200	-110	-216
250	-111	-217
300	-112	-218
400	-113	-219
500	-114	-220
600	-115	-221

