


CONNECTING TUBES

QC-6090 CONNECTING TUBE POLYPROPYLENE

"Y" shape, similar to QC-6085, but in a Y shape. Packs of 12.

Tube, O.D.


MM	Code
8	-110
9.5	-111
12.7	-112



QC-6095 CONNECTING TUBE POLYETHYLENE

Doubles as a stopcock and quick disconnect. A 90° twist turns it on and off. Approximately 83mm long, fits tubing 6mm-8mm ID.

Description	Code
Pack (6)	-110
Case (36)	-111




CONDUCTIVITY CELLS

QC-7000 CONDUCTIVITY CELL

Low, borosilicate glass, one arm on each side of 30ml center well, 100mm high. Complete with platinum electrodes, thermometer and stopper.

Code
-110



QC-7020 CONDUCTIVITY CELL

Medium, borosilicate glass, one arm on each side, 35ml, 80mm high. Complete with platinum electrodes, thermometer and stopper.

Code
-110

