

CENTRIFUGE TUBES


QC-6010 CENTRIFUGE TUBE

Pear shaped, 100ml with a 1.5ml stem.

Graduations, ML

0 to 1.5 x .1
 1.5 to 5 x .5
 5 to 10 x 1.0
 10 to 25 x 5
 and at 50 and 100

Code
 -110




QC-6020 CENTRIFUGE TUBE

Pear shaped, 100ml with a 1ml stem. Used for phosphorous determinations in iron and steel.

Graduations, ML

0 to 1 x 0.05
 1 to 5 x 1.0
 5 to 50 x 5
 50 to 100 x 10

Code
 -110



CONNECTING TUBES

QC-6050 CONNECTING TUBE

Straight, hose connections on both ends for rubber or plastic hose fittings.


| Tube O.D. MM | Middle of Hose Connection O.D., MM | Code |
|--------------|------------------------------------|------|
| 8 | 8.5 | -110 |
| 10 | 11 | -111 |
| 12 | 14 | -112 |



QC-6080 CONNECTING TUBE

"Y" shape. Similar to QC-6060 but in a Y shape.


| Tube O.D. MM | Middle of Hose Connection O.D., MM | Code |
|--------------|------------------------------------|------|
| 8 | 8.5 | -110 |
| 10 | 11 | -111 |
| 12 | 14 | -112 |



QC-6060 CONNECTING TUBE

"T" shape. Similar to QC-6050 but with a third connection on side.


| Tube O.D. MM | Middle of Hose Connection MM | Code |
|--------------|------------------------------|------|
| 8 | 8.5 | -110 |
| 10 | 11 | -111 |
| 12 | 14 | -112 |



QC-6085 CONNECTING TUBE POLYPROPYLENE

"T" shape with hose connections on all ends for rubber or plastic fittings. Autoclavable and has good chemical resistance. Packs of 12.

| Tube, O.D., MM | Code |
|----------------|------|
| 8 | -110 |
| 9.5 | -111 |
| 12.7 | -112 |




CONNECTING TUBES

QC-6090 CONNECTING TUBE POLYPROPYLENE

"Y" shape, similar to QC-6085, but in a Y shape. Packs of 12.

Tube, O.D.


| MM | Code |
|------|------|
| 8 | -110 |
| 9.5 | -111 |
| 12.7 | -112 |



QC-6095 CONNECTING TUBE POLYETHYLENE

Doubles as a stopcock and quick disconnect. A 90° twist turns it on and off. Approximately 83mm long, fits tubing 6mm-8mm ID.

| Description | Code |
|-------------|------|
| Pack (6) | -110 |
| Case (36) | -111 |




CONDUCTIVITY CELLS

QC-7000 CONDUCTIVITY CELL

Low, borosilicate glass, one arm on each side of 30ml center well, 100mm high. Complete with platinum electrodes, thermometer and stopper.

| Code |
|------|
| -110 |



QC-7020 CONDUCTIVITY CELL

Medium, borosilicate glass, one arm on each side, 35ml, 80mm high. Complete with platinum electrodes, thermometer and stopper.

| Code |
|------|
| -110 |

