


GAS APPARATUS

QG-2115 GAS BALLOON

With one 0-4 range high vacuum teflon threaded valve.


Capacity, ML	Code
100	-110
250	-111
500	-112
1000	-113



QG-2120 GAS BULB

With break seal at the top and filling tube on bottom constricted for easy seal off.

Capacity ML	Tube OD, MM	Code
50	10	-110
100	12	-111
250	12	-112
500	15	-113
1000	15	-114



REMEMBER

"Q" Glass has been working with customers to design and build custom glass apparatus and component parts since 1946. We can do the same for you. Please write or call us with your special requirements.

GRADED SEALS

QG-2140 GRADED SEAL

Kovar to borosilicate.

Kovar Diameter, Inches	Code
1/8	-110
1/4	-111
3/8	-112
1/2	-113
3/4	-114
1	-115

QG-2150 GRADED SEAL

Stainless steel to borosilicate.

Stainless steel Diameter, Inches	Code
1/8	-110
1/4	-111
3/8	-112
1/2	-113
3/4	-114
1	-115



QG-2160 GRADED SEAL

Quartz to borosilicate. Each grade has approximately 3 inches of quartz and 4 inches of borosilicate.

Quartz ID x OD	Code	Quartz ID x OD	Code
2 x 4	-110	16 x 18	-223
3 x 5	-111	17 x 19	-224
4 x 6	-112	18 x 20	-225
5 x 7	-113	20 x 22	-226
6 x 8	-114	22 x 25	-227
7 x 9	-115	25 x 28	-228
8 x 10	-116	27 x 30	-229
10 x 12	-117	30 x 33	-230
11 x 13	-118	32 x 35	-231
12 x 14	-119	35 x 38	-232
13 x 15	-120	37 x 40	-233
14 x 16	-121	40 x 43	-234
15 x 17	-122		

