

HYDROMETERS

QH-5070 HYDROMETER FOR PETROLEUM

Hydrometer, API Scale, 1/10° Divisions,
335mm Length

Range	Code	Range	Code
-1 to 11°	-110	49 to 61°	-230
9 to 21°	-111	59 to 71°	-231
19 to 31°	-112	69 to 81°	-232
29 to 41°	-113	79 to 91°	-233
39 to 51°	-114	89 to 101°	-234

Hydrometer, API Scale, 1/10° Divisions,
169mm length

Range	Code	Range	Code
0 to 6°	-115	50 to 56°	-235
5 to 11°	-116	55 to 61°	-236
10 to 16°	-117	60 to 66°	-237
15 to 21°	-118	65 to 71°	-238
20 to 26°	-119	70 to 76°	-239
25 to 31°	-120	75 to 81°	-240
30 to 36°	-121	80 to 86°	-241
35 to 41°	-122	85 to 91°	-242
40 to 46°	-123	90 to 96°	-243
45 to 51°	-124	95 to 101°	-244

Hydrometer, API Scale, 2/10 Divisions,
170mm length

Range	Code	Range	Code
-1 to 11°	-125	49 to 61°	-245
9 to 21°	-126	59 to 71°	-246
19 to 31°	-127	69 to 81°	-247
29 to 41°	-128	79 to 91°	-248
39 to 51°	-129	89 to 101°	-249



QH-5075 THERMOHYDROMETER FOR PETROLEUM

Thermohydrometer, API Scale,
1/10° Divisions, 385mm Length.
Thermometer in body range
0 to 150°F, 2°F Divisions

Range	Code
15 to 23°	-110
22 to 30°	-111
29 to 37°	-112
36 to 44°	-113
43 to 51°	-114

Thermohydrometer, API Scale,
1/10° Divisions, 385mm length.
Thermometer in body range
0 to 150°F, 2°F Divisions

Range	Code	Range	Code
-1 to 11°	-115	49 to 61°	-120
9 to 21°	-116	59 to 71°	-121
19 to 31°	-117	69 to 81°	-122
29 to 41°	-118	79 to 91°	-123
39 to 51°	-119	89 to 101°	-124

Thermohydrometer, API Scale,
1/10° Divisions, 385mm length.
Thermometer in body range
30 to 220°F, 2°F Divisions

Range	Code
-1 to 11°	-125
9 to 21°	-126
19 to 31°	-127
29 to 41°	-128

Thermohydrometer, API Scale,
2/10° Divisions, 190mm length.
Thermometer in body range
20 to 130°F, 2°F Divisions

Range	Code	Range	Code
-1 to 11°	-129	49 to 61°	-134
9 to 21°	-130	59 to 71°	-135
19 to 31°	-131	69 to 81°	-136
29 to 41°	-132	79 to 91°	-137
39 to 51°	-133	89 to 101°	-138



