

MICRO LAB • BULBS

QM-500 BULB, KJELDAHL

Outer joint at top with inner joint below a round trap.


Joint Size	Code
14/20	-10
19/22	-11



QM-510 BULB, KJELDAHL

Inner joints of the same size on the arm and bottom. Bottom inner tube has opening for return of trapped liquid.

Joint Size	Code
14/20	-10
19/22	-11




MICRO LAB • CONDENSERS

QM-520 CONDENSER, ALLIHN

The top of the condenser is a plain tube, the bottom has a drip tip.

Jacket	14/20	19/22
Length, MM	Code	Code
150	-10	-22
200	-11	-23



QM-550 CONDENSER, AIR

Plain tube with drip tip at bottom. Tube is approximately 200mm long.


Joint Size	Code
14/20	-10
19/22	-11



QM-530 CONDENSER, ALLIHN

Similar to QM-520, but with an outer joint at top.


Jacket	14/20	19/22
Length, MM	Code	Code
150	-10	-22
200	-11	-23



QM-560 CONDENSER, LIEBIG

Top outer and bottom inner joints are the same size.


Jacket	14/20	19/22
Length, MM	Code	Code
110	-10	-22
250	-11	-23



QM-540 CONDENSER, ALLIHN

Full length, water cooled, inner and outer joints. The jacket is approximately 150mm.

Joint Size	Code
14/20	-10
19/22	-11



QM-570 CONDENSER, WEST

Top outer and bottom inner joints are the same size.

Jacket	14/20	19/22
Length, MM	Code	Code
110	-10	-22
200	-11	-23

