


MICRO LAB • DISTILLING HEADS

QM-780 DISTILLING HEAD


Removable cold finger is approximately 100mm long with a 14/20 inner joint. The 10/18 outer joint on main tube is for thermometer. Inner joints are the same size. Available with a glass or teflon stopcock.



Description	Inner Joints	Glass Stopcock Code	Teflon Stopcock Code
Complete	14/20	-10	-24
Complete	19/22	-11	-25
Head Only	14/20	-12	-26
Head Only	19/22	-13	-27
Cold Finger Only	14/20	QM-590-11	

QM-790 DISTILLING HEAD

Built in Allihn condenser and stopcock manifold for control of reflux and vacuum, with vent next to the oblique stopcock. 10/18 outer joint for thermometer and two inner joints of the same size. Available with glass or teflon stopcocks.




Inner Joints	Glass Stopcock Code	Teflon Stopcock Code
14/20	-10	-22
19/22	-11	-23

MICRO LAB • DISTILLING RECEIVERS

QM-800 DISTILLING RECEIVER


Graduated, conical shape bottom with outer joint and hooks. Graduated from 0-2ml in 0.1ml and 2-12 (15)ml in 0.2ml subdivisions.



Capacity, ML	14/20 Code	19/22 Code
12	-10	-22
15	-11	-23

QM-820 DISTILLING RECEIVER

Receiver has an outlet tube for vacuum with a stopcock and a 10ml graduated section with stopcock and drip tip at bottom. Available with glass or teflon stopcocks.




Joint Size	Glass Stopcock Code	Teflon Stopcock Code
14/20	-10	-22
19/22	-11	-23

QM-810 DISTILLING RECEIVER

With outer joint above the graduated section and inner joint of the same size on the side arm. 2ml has a cylindrical shaped reservoir tube and 5 and 10ml have conical shaped tubes. Available with glass or teflon stopcocks.

Capacity, ML	Subdivisions, ML	Joint Size	Glass Stopcock Code	Teflon Stopcock Code
2	0.02	14/20	-10	-26
5	0.2	14/20	-11	-27
10	0.1	14/20	-12	-28
2	0.02	19/22	-13	-29
5	0.2	19/22	-14	-30
10	0.1	19/22	-15	-31




MICRO LAB • DISTILLING RECEIVERS

QM-830 DISTILLING RECEIVER

Distributing chamber with four 15ml round bottom flasks as receivers. Top vacuum adapter is on a 105° angle with a long delivery tip. All joints are the same size.


Description	14/20 Code	19/22 Code
Complete	-10	-24
Distribution chamber only	-11	-25
Vacuum adapter only	-12	-26
15ml Flask only	-13	-27



QM-840 DISTILLING RECEIVER

Similar to QM-830, but the vacuum adapter is straight, with a 2mm teflon stopcock. All joints are the same size.

Description	14/20 Code	19/22 Code
Complete	-10	-24
Distribution chamber only	-11	-25
Vacuum adapter only	-12	-26
15ml Flask only	-13	-27



MICRO LAB • EXTRACTION APPARATUS

QM-850 EXTRACTION APPARATUS For use with a 10 x 50 glass or paper extraction thimble. Consists of an extractor body, a condenser and a 25ml round bottom flask. Joints between extractor and condenser are 19/22 and joints between extractor and flask are 14/20.

Description	Code
Complete	-10
Extractor Only	-11
Condenser Only	-12
25ml Flask Only	-13




MICRO LAB • FLASKS

QM-860 FLASK, SINGLE NECK

Round bottom with outer joint. Short neck.

Capacity	14/20	19/22	Capacity	14/20	19/22
ML	Code	Code	ML	Code	Code
2	-10	-25	50	-30	-44
5	-11	-26	100	-31	-45
10	-12	-27	250	-32	-46
15	-13	-28	500	-33	-47
25	-14	-29			



QM-870 FLASK, SINGLE NECK

Round bottom with outer joint. Long neck.

Capacity	14/20	19/22
ML	Code	Code
5	-10	-25
10	-11	-26
25	-12	-27
50	-13	-28
100	-14	-29

