


# MICRO LAB • DISTILLING RECEIVERS

## QM-830 DISTILLING RECEIVER

Distributing chamber with four 15ml round bottom flasks as receivers. Top vacuum adapter is on a 105° angle with a long delivery tip. All joints are the same size.


| Description               | 14/20<br>Code | 19/22<br>Code |
|---------------------------|---------------|---------------|
| Complete                  | -10           | -24           |
| Distribution chamber only | -11           | -25           |
| Vacuum adapter only       | -12           | -26           |
| 15ml Flask only           | -13           | -27           |



## QM-840 DISTILLING RECEIVER

Similar to QM-830, but the vacuum adapter is straight, with a 2mm teflon stopcock. All joints are the same size.

| Description               | 14/20<br>Code | 19/22<br>Code |
|---------------------------|---------------|---------------|
| Complete                  | -10           | -24           |
| Distribution chamber only | -11           | -25           |
| Vacuum adapter only       | -12           | -26           |
| 15ml Flask only           | -13           | -27           |



# MICRO LAB • EXTRACTION APPARATUS

**QM-850 EXTRACTION APPARATUS** For use with a 10 x 50 glass or paper extraction thimble. Consists of an extractor body, a condenser and a 25ml round bottom flask. Joints between extractor and condenser are 19/22 and joints between extractor and flask are 14/20.

| Description     | Code |
|-----------------|------|
| Complete        | -10  |
| Extractor Only  | -11  |
| Condenser Only  | -12  |
| 25ml Flask Only | -13  |




# MICRO LAB • FLASKS

## QM-860 FLASK, SINGLE NECK

Round bottom with outer joint. Short neck.

| Capacity | 14/20 | 19/22 | Capacity | 14/20 | 19/22 |
|----------|-------|-------|----------|-------|-------|
| ML       | Code  | Code  | ML       | Code  | Code  |
| 2        | -10   | -25   | 50       | -30   | -44   |
| 5        | -11   | -26   | 100      | -31   | -45   |
| 10       | -12   | -27   | 250      | -32   | -46   |
| 15       | -13   | -28   | 500      | -33   | -47   |
| 25       | -14   | -29   |          |       |       |



## QM-870 FLASK, SINGLE NECK

Round bottom with outer joint. Long neck.

| Capacity | 14/20 | 19/22 |
|----------|-------|-------|
| ML       | Code  | Code  |
| 5        | -10   | -25   |
| 10       | -11   | -26   |
| 25       | -12   | -27   |
| 50       | -13   | -28   |
| 100      | -14   | -29   |

