

MICRO LAB • FUNNELS

QM-1065 FUNNEL, SEPARATORY CYLINDRICAL

Similar to QM-1060, but with 2mm teflon stopcocks.

Capacity	Subdivisions	14/20	19/22
ML	ML	Code	Code
10	0.2	-10	-24
25	0.5	-11	-25
50	1.0	-12	-26
60	1.0	-13	-27



QM-1080 FUNNEL, SEPARATORY, SQUIBB

Pear shaped funnel with an outer joint at top and an inner joint with drip tube at bottom below a 2mm glass stopcock. Complete with stopper.

Capacity, ML	14/20 Code	19/22 Code
10	-10	-25
25	-11	-26
30	-12	-27
60	-13	-28
125	-14	-29



QM-1070 FUNNEL, SEPARATORY, SQUIBB

Pear shaped funnel with an outer joint at top and a 2mm glass stopcock at bottom. Complete with stopper.

Capacity	14/20	19/22
ML	Code	Code
10	-10	-25
25	-11	-26
30	-12	-27
60	-13	-28
125	-14	-29



QM-1085 FUNNEL, SEPARATORY, SQUIBB

Similar to QM-1080, but with a 2mm teflon stopcock.

Capacity	14/20	19/22
ML	Code	Code
10	-10	-25
25	-11	-26
30	-12	-27
60	-13	-28
125	-14	-29



QM-1075 FUNNEL, SEPARATORY, SQUIBB

Similar to QM-1070, but with a 2mm teflon stopcock.


Capacity	14/20	19/22
ML	Code	Code
10	-10	-25
25	-11	-26
30	-12	-27
60	-13	-28
125	-14	-29



QM-1090 FUNNEL, POWDER

Short form funnel with an inner joint at bottom. Top diameter is 50mm.

Joint Size	Code
14/20	-10
19/22	-11




MICRO LAB • STIRRERS

QM-2000 STIRRER BEARING 6MM

6mm ID precision bore, With mercury cup at top and inner joint at bottom.

Joint Size	Code
14/20	-10
19/22	-11




MICRO LAB • STIRRERS

QM-2010 STIRRER BEARING 10MM

10mm ID precision bore, with mercury cup at top and inner joint at bottom.

Joint Size	Code
14/20	-10
19/22	-11



QM-2040 STIRRER ROD 6MM

6mm OD x 300mm long. The entire length is precision ground and has bottom button for attachment of stirrer blades.


Code
-10



QM-2020 STIRRER BEARING JACKETED 6MM

6mm ID precision bore, jacketed, with inlet and outlet hose connections.

Joint Size	Code
14/20	-10
19/22	-11



QM-2050 STIRRER ROD 10MM

10mm OD x 400mm long. The entire length is precision ground and has bottom button for attachment of stirrer blades.


Code
-10



QM-2030 STIRRER BEARING JACKETED 10MM

10mm ID precision bore, jacketed, with inlet and outlet hose connections.

Joint Size	Code
14/20	-10
19/22	-11



QM-2060 STIRRER BLADES, TEFLON

Crescent type stirrer blade made of solid teflon. To fit 6mm or 10mm stirrer rods.

Stirrer, OD, MM	Blade Length, MM	Code
6	33	-10
10	46	-11



MICRO LAB • STOPPERS

QM-2070 STOPPER, PENNYHEAD

Hollow stoppers with a closed bottom, ground for use with outer joints.

Joint Size	Code
10/18	-10
14/20	-11
19/22	-12

