


POUR FOAM TEST APPARATUS

QP-3090 POUR FOAM TEST APPARATUS

Ross-Miles, ASTM #D1173, used for testing the foaming properties of surface-active agents. Consists of a 200ml capacity foam pipet with 2mm teflon stopcock and a foam receiver graduated at 50 and 250ml with a 6mm teflon stopcock.

<u>Description</u>	<u>Code</u>
Complete	-110
Foam Pipet Only	-111
Foam Receiver	-112



PROCESS PIPE

QP-4000 PROCESS PIPE

Heavy wall glass pipe. Available conical or beaded fittings.

All sizes available in straight lengths, elbows, reducers, T-tubes, hose connection adapters, thermometer wells, etc . . . Clamps and gaskets to fit are also available.

Please call with your requirements.

