

ULTRASONIC CLEANERS

QU-1000 ULTRASONIC CLEANER

Industrial cleaning performance and reliability in a benchtop cleaner. Crest's exclusive tru-sweep feature ensures uniform cleaning throughout the tank by sweeping the ultrasonic frequency + 2kHz, creating overlapping ultrasonic waves. This eliminates inconsistent cleaning. Timer has continuous operation capabilities. Metal skirts resist damage from chemical spills and will not crack. Drains with flow control valve on 3/4 gallon units and larger. Choose from two models, one with a 0-30 minute timer or one with a 0-30 minute timer and an adjustable, thermostatically controlled heater (ambient 80° C). Complete with cover on all units



Tank Capacity Gallons/Liter	Inside Tank Dimensions L x W x D In.	Overall Dimensions L x W x D In.	Sonic Power Avg/Peak	Timer Code	Heat & Timer Code
1/2 1.75	5 3/4 x 5 1/4 x 4	6 1/4 x 5 3/4 x 9	45 135	-110	-217
3/4 2.75	9 1/4 x 5 1/4 x 4	9 3/4 x 5 3/4 x 9	90 270	-111	-218
1 1/2 5.75	11 1/2 x 6 x 6	12 1/8 x 6 1/4 x 11	135 405	-112	-219
1 3/4 6.9	19 1/2 x 5 1/2 x 6	20 3/8 x 5 3/4 x 11	180 540	-113	-220
2 1/2 9.5	11 1/2 x 9 1/2 x 6	12 1/8 x 5 3/4 x 11	180 540	-114	-221
5 1/2 18.75	19 1/2 x 11 1/2 x 6	20 3/8 x 12 1/8 x 13	270 810	-115	-222
7 1/2 28.0	19 1/2 x 11 1/2 x 8	20 3/8 x 12 1/8 x 15	360 1080	-116	-223

QU-2000 ULTRASONIC CLEANER, DIGITAL

Similar to the dial units, but completely digital with precise control of time (0-99 minutes), heat (ambient 80° C). power level and degas functions. Complete with cover on all sizes.



Tank Capacity Gallons/Liter	Inside Tank Dimensions L x W x D In.	Overall Dimensions L x W x D In.	Sonic Power Avg/Peak	Code
1/2 1.75	5 3/4 x 5 1/4 x 4	6 1/4 x 5 3/4 x 9	45 135	-110
3/4 2.75	9 1/4 x 5 1/4 x 4	9 3/4 x 5 3/4 x 9	90 270	-111
1 1/2 5.75	11 1/2 x 6 x 6	12 1/8 x 6 1/4 x 11	135 405	-112
1 3/4 6.9	19 1/2 x 5 1/2 x 6	20 3/8 x 5 3/4 x 11	180 540	-113
2 1/2 9.5	11 1/2 x 9 1/2 x 6	12 1/8 x 5 3/4 x 11	180 540	-114
5 1/2 18.75	19 1/2 x 11 1/2 x 6	20 3/8 x 12 1/8 x 13	270 810	-115
7 1/2 28.0	19 1/2 x 11 1/2 x 8	20 3/8 x 12 1/8 x 15	360 1080	-116

QU-3000 ULTRASONIC CLEANING SOLUTIONS

Choose from three cleaning solutions. All non-ammoniated phosphate free and biodegradable. Available in one gallon containers or cases of 4 one gallon containers.

Description	Each Code	Case Code
• General purpose cleaner, recommended for labware and jewelry -10:1 dilution with water.	-110	-213
• Oxide remover (Copper and Brass)	-111	-214
• Buffing compound remover (not recommended for aluminum)	-112	-215

Special requirements for all of your ultrasonic needs, such as custom sizes, options, etc. are available. Please call or fax us with your requirements.