

# BEAKERS

## GRIFFEN LOW FORM GRADUATED BEAKER WITH POUR SPOUT

<u>Part No.</u>	<u>Capacity</u>	<u>Pack</u>	<u>Price</u>
QBS-100-10	50ml	12/pack	\$35.25
QBS-100-11	100ml	12/pack	\$37.75
QBS-100-12	150ml	12/pack	\$37.75
QBS-100-13	250ml	12/pack	\$45.25
QBS-100-14	500ml	12/pack	\$27.35
QBS-100-18	600ml	6/pack	\$34.25
QBS-100-15	800ml	6/pack	\$29.75
QBS-100-16	1000ml	6/pack	\$51.65
QBS-100-17	2000ml	4/pack	\$77.80



# BOTTLES

## REAGENT BOTTLE, OUTER JOINT AT TOP COMPLETE WITH GLASS STOPPER


<u>Part No.</u>	<u>Capacity</u>	<u>Case</u>	<u>Price</u>
QBS-105-10	125ml	12/case	\$82.05
QBS-105-11	250ml	12/case	\$94.20
QBS-105-12	500ml	12/case	\$108.20
QBS-105-13	1000ml	12/case	\$103.35



# LABORATORY BURETS

## LABORATORY GRADE GRADUATED BURET 2MM TEFLON STOPCOCK


<u>Part No.</u>	<u>Capacity</u>	<u>Subdivisions</u>	<u>Tolerance</u>	<u>Price</u>
QBS-126-10	10ml	0.05ml	0.05ml	\$43.15
QBS-126-11	25ml	0.1ml	0.1ml	\$49.85
QBS-126-12	50ml	0.1ml	0.1ml	\$57.75
QBS-126-13	100ml	0.1ml	0.1ml	\$66.25



# CLASS A BURETS

## CLASS A, GRADUATED BURET 2MM TEFLON STOPCOCK


<u>Part No.</u>	<u>Capacity</u>	<u>Subdivisions</u>	<u>Tolerance</u>	<u>Price</u>
QBS-127-10	10ml	0.05ml	0.02ml	\$89.35
QBS-127-11	25ml	0.1ml	0.03ml	\$95.45
QBS-127-12	50ml	0.1ml	0.05ml	\$101.50
QBS-127-13	100ml	0.2ml	0.1ml	\$120.35



# LABORATORY CYLINDERS

**GRADUATED, LABORATORY GRADE,  
CALIBRATED TO DELIVER, WITH POUR SPOUT**


<u>Part No.</u>	<u>Capacity</u>	<u>Subdivision</u>	<u>Pack</u>	<u>Price</u>
QCS-703-10	25ml	0.5ml	4/pack	\$29.75
QCS-703-11	50ml	1ml	4/pack	\$36.40
QCS-703-12	100ml	1ml	4/pack	\$43.15
QCS-703-13	250ml	2ml	2/pack	\$36.40
QCS-703-14	500ml	5ml	1/pack	\$27.35
QCS-703-15	1000ml	10ml	1/pack	\$36.40
QCS-703-16	2000ml	20ml	1/pack	\$50.20



# CLASS A CYLINDERS

**GRADUATED, CLASS A,  
CALIBRATED TO DELIVER, WITH POUR SPOUT**

<u>Part No.</u>	<u>Capacity</u>	<u>Subdivisions</u>	<u>Tolerance</u>	<u>Pack</u>	<u>Price</u>
QCS-704-10	50ml	1.0ml	0.25ml	4/pack	\$21.20
QCS-704-11	100ml	1.0ml	0.5ml	4/pack	\$26.50
QCS-704-12	250ml	2.0ml	1ml	4/pack	\$31.85
QCS-704-13	500ml	5.0ml	2ml	2/pack	\$48.55
QCS-704-14	1000ml	10.0ml	3ml	1/pack	\$66.85



# CONDENSERS

## ALLIHN BULB TYPE


300MM EFF. LENGTH  
24/40 JOINTS

Part No.

QCS-300-10

Price

\$ 63.20



## GRAHAM COIL TYPE


300MM EFF. LENGTH  
24/40 JOINTS.

Part No.

QCS-308-10

Price

\$ 79.00



## LIEBIG


300MM EFF. LENGTH  
24/40 JOINTS

Part No.

QCS-402-10

Price

\$ 59.55



## SPIRAL REFLUX


300MM EFF. LENGTH  
24/40 JOINTS.

Part No.

QCS-503-10

Price


\$ 108.20



# ERLENMEYER FLASKS

## GRADUATED


<u>Part No.</u>	<u>Capacity</u>	<u>Pack</u>	<u>Price</u>
QFS-307-10	100ml	12/pack	\$39.50
QFS-307-11	250ml	12/pack	\$42.50
QFS-307-12	500ml	6/pack	\$27.25
QFS-307-13	1000ml	6/pack	\$49.95
QFS-307-14	2000ml	1/pack	\$24.25



# VOLUMETRIC FLASKS

## CLASS A, COMPLETE WITH STOPPER


<u>Part No.</u>	<u>Capacity</u>	<u>Tolerance</u>	<u>Pack</u>	<u>Price</u>
QFS-406-10	50ml	0.05	6/pack	\$139.80
QFS-406-11	100ml	0.08	6/pack	\$139.80
QFS-406-12	200ml	0.10	6/pack	\$179.00
QFS-406-13	250ml	0.12	6/pack	\$203.50
QFS-406-14	500ml	0.20	6/pack	\$253.50
QFS-406-15	1000ml	0.30	1/pack	\$54.70
QFS-406-16	2000ml	0.50	1/pack	\$66.85



# SEPARATORY FUNNELS

## CYLINDRICAL, UNGRADUATED GLASS STOPCOCK, COMPLETE WITH STOPPER


<u>Part No.</u>	<u>Capacity</u>	<u>Price</u>
QFS-500-10	50ml	\$42.55
QFS-500-11	100ml	\$48.65
QFS-500-12	250ml	\$60.75
QFS-500-13	500ml	\$79.00
QFS-500-14	1000ml	\$115.50



# SEPARATORY FUNNELS

## CYLINDRICAL, GRADUATED GLASS STOPCOCK AT BOTTOM COMPLETE WITH STOPPER

<u>Part No.</u>	<u>Capacity</u>	<u>Price</u>
QFS-508-10	50ml	\$79.00
QFS-508-11	125ml	\$85.10
QFS-508-12	250ml	\$97.25
QFS-508-13	500ml	\$120.35
QFS-508-14	1000ml	\$164.10



# SEPARATORY FUNNELS

PEAR-SHAPED, STOPPER JOINT AT TOP  
TEFLON STOPCOCK AT BOTTOM


<u>Part No.</u>	<u>Capacity</u>	<u>Price</u>
QFS-604-10	125ml	\$75.30
QFS-604-11	250ml	\$83.90
QFS-604-12	500ml	\$89.95
QFS-604-13	1000ml	\$107.20
QFS-604-14	2000ml	\$114.25



# POWDER FUNNELS

SHORT STEM, 60°

<u>Part No.</u>	<u>Top Dia.</u>	<u>Stem Dia. x Length</u>	<u>Pack</u>	<u>Price</u>
QFS-701-10	75mm	9mm x 75mm	12/pack	\$59.55
QFS-701-11	100mm	13mm x 100mm	6/pack	\$65.65
QFS-701-12	125mm	15mm x 120mm	2/pack	\$43.75
QFS-701-13	150mm	17mm x 150mm	2/pack	\$90.55



# VOLUMETRIC PIPETS

## CLASS A, CALIBRATED TO DELIVER, COLOR CODED

<u>Part No.</u>	<u>Capacity</u>	<u>Tolerance</u>	<u>Color Code</u>	<u>Pack</u>	<u>Price</u>
QPS-302-10	5ml	0.015	white	12/pack	\$57.75
QPS-302-11	10ml	0.02	red	12/pack	\$60.20
QPS-302-12	20ml	0.03	yellow	1/pack	\$7.90
QPS-302-13	25ml	0.03	blue	1/pack	\$8.20
QPS-302-14	50ml	0.05	red	1/pack	\$10.65
QPS-302-15	100ml	0.08	yellow	1/pack	\$13.30



# MEASURING PIPETS

## CLASS A, CALIBRATED TO DELIVER, COLOR CODED

<u>Part No.</u>	<u>Capacity</u>	<u>Tolerance</u>	<u>Color Code</u>	<u>Pack</u>	<u>Price</u>
QPS-303-10	1ml	.01 / .006	yellow	6/pack	\$47.40
QPS-303-11	2ml	.02 / .01	black	6/pack	\$47.40
QPS-303-12	5ml	.05 / .03	red	6/pack	\$47.40
QPS-303-13	10ml	0.1 / .05	orange	6/pack	\$59.55
QPS-303-14	25ml	0.2 / .1	white	6/pack	\$59.55





# CENTRIFUGE TUBES

**BOROSILICATE GLASS, GRADUATED, CONICAL BOTTOM, BEADED TOP**

<u>Part No.</u>	<u>Capacity</u>	<u>Subdivisions</u>	<u>Price</u>
QCS-504-10	10ml	0.2ml	\$12.15
QCS-504-11	25ml	0.5ml	\$14.00
QCS-504-12	50ml	1ml	\$14.60



# TEST TUBES

**ROUND BOTTOM WITHOUT RIM**

<u>Part No.</u>	<u>O.D. x Length</u>	<u>Pack</u>	<u>Price</u>
QTS-304-10	10mm x 75mm	1,000/pack	\$54.70
QTS-304-11	12mm x 75mm	1,000/pack	\$72.95
QTS-304-12	13mm x 100mm	1,000/pack	\$85.10
QTS-304-13	16mm x 125mm	1,000/pack	\$115.50
QTS-304-14	16mm x 150mm	1,000/pack	\$127.65
QTS-304-15	18mm x 150mm	600/pack	\$102.10
QTS-304-16	25mm x 150mm	100/pack	\$97.25

